



Lympha Press USA

LYMPHA PRESS® 201MAX

HCPCS code E0652

- Single phase calibrated gradient sequential compression
- Quick 26-second cycle time
- Tamper-proof lock and automatic shutoff
- Unilateral or bilateral operation
- Use with Lympha Press® quality garments



Calibrated gradient sequential compression therapy system

LYMPHA PRESS® 201MAX is a home care system for treatment of lymphedema and venous insufficiency. It applies a gentle but effective sequential compression wave to maximally stimulate the lymphatic and venous systems as it moves fluids to healthy areas for uptake.

201MAX is easy to use and is fully adjustable, to meet individual patient needs. Pressure can be set according to physician prescription: low pressure for light massage, or higher pressure for those whose condition requires it. Its three zone gradient allows for more intense therapy where needed, while minimizing pressure over areas of discomfort. In addition, our narrow, overlapping air chambers can be individually disconnected without disrupting therapy to adjacent areas.



Indications:

- Primary lymphedema
- Secondary lymphedema
- Venous insufficiency
- Venous stasis ulcers
- Dysfunction of the muscle pump

Contraindications:

- Known or suspected DVT or pulmonary embolism
- Inflammatory phlebitis
- Acute infection
- Decompensated cardiac failure
- Severe arteriosclerosis or other ischemic vascular disease
- Any circumstance in which increased venous and lymphatic return is undesirable

Setting the pressure on the 201MAX is quick, easy and accurate, and can be calibrated for every zone, size and garment individually. The tamper-proof mechanism and auto off switch increase patient compliance



Features:

- Effective single-phase calibrated gradient compression, with 10-24 overlapping chambers and three individual zones of pressure. Our pressure mechanism prevents inadvertent setting of a reverse gradient.
- Lympha Press® quick cycle (26 seconds) for maximum therapy and patient comfort.
- Easy to read digital pressure settings can be fine-tuned for every patient individually.
- Tamper-proof knobs and automatic treatment time for patient compliance.
- Unilateral or bilateral operation.



Technical Specifications:

Model	201Max Sequential
Cycle	26 sec inflate/4 sec deflate
Pressure Range	20 – 80 mmHg
Size	6 × 12.5 × 15.3 inches
Weight	11.7 lb.
Electrical	115V, 60 Hz

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Lympha Press USA

Distributed by: Lympha Press USA Ltd. 232 Park Avenue, Manalapan, NJ 07726
 Tel: 1-888-596-7421, 732-792-9677, Fax: 732-792-9745
 info@lymphapress.com www.lymphapress.com