


# Lympha Press USA

Website: [www.lymphapress.com](http://www.lymphapress.com)

232 Park Avenue Manalapan New Jersey, 07726  
 1-888-596-7421 ♦ 732-792-9677 ♦ 732-792-9710 ♦ 732-792-9743  
 Fax: 732-792-9745

## LymphaPod® Garment



Dimensions are in cm.

LymphaPod® Left Hose Exit Cat. No. L4MP800L LymphaPod® Right Hose Exit Cat. No. L4MP800R					
Reference Point	On patient (approximate*)	Circumference	With one expander Type A (one side only**)	With two expanders Type A (one on each side)	Addition per expander Type A
A	Waist	206	222	238	16
B	Hips	212	227	242	15
C	Upper Thigh	130	143	143	13
D	Knee	125	138	138	13
E	Ankle/Heel	104	114	114	10
F	End	87	95	95	8

\* Provided for illustrative purposes only. Location of patient landmarks will vary depending on height and leg length. Contact your Lympha Press® dealer for assistance.

\*\* Leg measurements for side with expander added.