


# Lympha Press USA

Website: [www.lympha-press.com](http://www.lympha-press.com)

232 Park Avenue Manalapan New Jersey, 07726  
 1-888-596-7421 • 732-792-9677 • 732-792-9710 • 732-792-9743  
 Fax: 732-792-9745

## Comfy Sleeves and Expanders



Comfy sleeve.vsd

### Notes:

1. All Dimensions are in cm (Centimeters).
2. General tolerance:  $\pm 1$ cm
3. A, C & D dimensions are circumferences.

Cat #	Sleeve details and length			Sleeve's dimensions					Expander Cat. No.	Expander sizes	Dimensions with Expander		
				A	B	C	D	L			A	C	D
L41L2285	Comfy	M	85cm	70	30	32	40	83	L415A087	Size C	80	42	50
L41L2365	Comfy	L	65cm	80	40	47	55	65	L415A078	Size B	90	57	65
L41L2375	Comfy	L+	75cm	90	40	47	55	75	L415A087	Size C	100	57	65
L41L2385	Comfy	L++	85cm	94	40	47	55	84	L415A094	Size D	100	57	65