


# Lympha Press USA


Website: [www.lympha-press.com](http://www.lympha-press.com)

232 Park Avenue Manalapan New Jersey, 07726  
 1-888-596-7421 • 732-792-9677 • 732-792-9710 • 732-792-9743  
 Fax: 732-792-9745

## Lympha Pants II Dimensional data



## Expanders and Adapter



General tolerance:  $\pm 2$  cm

Garment Measuring Point	Bilateral pants				
	No Expanders	1 expander one side only	2 expanders one per side	4 expanders two per side	Increase per additional expander
Waist (A)	130.0	148.0	166.0	202.0	18.0
Hips at crotch (B)	127.0	142.0	157.0	187.0	15.0
Upper Thigh (C & D)	73.5	87.5	87.5	101.5	14.0
E (1/3 from end)	63.0	75.0	75.0	87.0	12.0
F (end)	52.5	61.5	61.5	70.5	9.0

Garment Measuring Point	Unilateral pants				
	No Expanders	1 expander one side only	2 expanders one per side	4 expanders two per side	Increase per additional expander
Waist (A)	124.5	142.5	160.5	196.5	18.0
Hips at crotch (B)	121.5	136.5	151.5	181.5	15.0
Upper Thigh, untreated leg (C)	68.0	68.0*	82.0**	96.0**	14.0
Upper Thigh, treated leg (D)	73.5	87.5	87.5	101.5	14.0
E (1/3 from end)	63.0	75.0	75.0	87.0	12.0
F (end)	52.5	61.5	61.5	70.5	9.0

\*assumes single Leg (long) expander is used on treated side

\*\* assumes GA (short) expander is used on un treated side