


lympho Press USA

Website: www.lympho-press.com

232 Park Avenue Manalapan New Jersey, 07726
 1-888-596-7421 • 732-792-9677 • 732-792-9710 • 732-792-9743
 Fax: 732-792-9745

PBS Sleeves with Tapered Expanders in cm



Notes:

1. All dimensions are in centimeters (cm)
2. General tolerance: ± 2 cm

Garment Type		Dimensions								Dimensions with Expander				
Cat No.	Size	A	B	C	D	E	F	L	Size	A	B	C	D	E
2250	Leg S	57	54	46	43	44	26	50	AT	67	63	50	48	48
2275	Leg S+	66	56	46	43	44	26	73	BT	76	66	54	48	48
2285	Leg M	70	60	47	43	44	26	84	CT	80	70	55	48	48
2295	Leg M+	80	72	60	48	37	32	95	ET	90	82	72	53	41
2365	Leg L	81	73	60	56	54	37	65	BT	91	83	70	59	58
2375	Leg L+	85	76	61	56	54	37	75	CT	95	86	69	59	58
2385	Leg L++	90	80	63	56	54	37	85	DT	100	90	71	59	58
1290	Arm M	58	--	--	29	--	--	88	--	--	--	--	--	--