


# Lympha Press USA

Website: [www.lympha-press.com](http://www.lympha-press.com)

232 Park Avenue Manalapan New Jersey, 07726  
 1-888-596-7421 • 732-792-9677 • 732-792-9710 • 732-792-9743  
 Fax: 732-792-9745

## PBS Sleeves and Expanders



1. All Dimensions are in cm (Centimeters).
2. General tolerance:  $\pm 2$  cm
3. A, C & D dimensions are circumferences.

Sleeve details and length	Sleeve's dimensions						Straight Expander				Tapered Expander							
							Sleeve's dimensions with Expander				Sleeve's dimensions with Expander							
	A	B	C	D	G	L	Cat. #	Size	A	C	D	G	Cat. #	Size	A	C	D	G
Leg S 50cm	56	25	32	40	43	52	L415A056	A	66	42	50	53	L415T056	AT	66	36	47	47
Leg S+ 75cm	64	25	32	40	43	72	L415A078	B	74	42	50	53	L415T078	BT	74	36	47	47
Leg M 85cm	68	25	32	40	43	81	L415A087	C	78	42	50	53	L415T087	CT	78	36	47	47
Leg M+ 90cm	72	25	32	40	43	82	L415A094	F	82	42	50	53	L415T094	FT	82	36	47	47
Leg L 65cm	80	33	47	53	52	70	L415A078	B	90	57	63	62	L415T078	BT	90	51	60	56
Leg L+ 75cm	86	33	47	53	52	79	L415A087	C	96	57	63	62	L415T087	CT	96	51	60	56
Leg L++ 85cm	88	33	47	53	52	87	L415A097	D	98	57	63	62	L415T097	DT	98	51	60	56
Arm M 90cm	56	--	29	--	--	88	L415A087	C	66	39	--	--	--	--	--	--	--	--