



Lympha Press USA

the Petite Basic System



HCPCS Code E0651



- Four fully overlapping compression cells
- Rapid cycle time
- Adjustable pause
- Lock feature to prevent tampering



The Petite Basic System™

Manufactured by the makers of Lympha Press®, **The Petite Basic System™** is a sequential pneumatic compression system specially designed for optimal reduction and management of peripheral edema, lymphedema, and ulceration associated with venous insufficiency.

The Petite Basic System™ has *four fully overlapping compression cells*

to ensure distal to proximal flow and comfortable compression. Pressure and pause times are easy to adjust. The *Lock feature* prevents accidental tampering, for patient safety.

The **Petite Basic System™** garments have a zipper closure for easy wear and are available in a wide range of sizes, with expanders for a temporary increase in circumference.



Indications

- Venous stasis ulcers
- Venous insufficiency
- Reduction and management of peripheral edema
- Lymphatic disorders

the PetiteBasicSystem



Garments Measurements in cm

Cat. No.	Garment Size & Length	Garment Dimensions								With Expander				
		A	B	C	D	E	F	L	Expander Size	A	B	C	D	E
2250	Leg S 50 cm	57	54	46	43	44	26	50	AT	67	63	50	48	48
2275	Leg S+ 73 cm	66	56	46	43	44	26	73	BT	76	66	54	48	48
2285	Leg M 84 cm	70	60	47	43	44	26	84	CT	80	70	55	48	48
2365	Leg L 65 cm	81	73	60	56	54	37	65	BT	91	83	70	59	58
2375	LegL+ 75 cm	85	76	61	56	54	37	75	CT	95	86	69	59	58
2385	Leg L++ 85 cm	90	80	63	56	54	37	85	DT	100	90	71	59	58
1290	Arm M 88 cm	58	-	-	29	-	-	88	C	68	-	-	39	-


A: Sleeve upper circumference, **L=** Sleeve length. **B - F** are additional measurements.

The Petite Basic System™

is designed to give quality treatment at home. Its ease of use encourages patient compliance, while its unique overlapping chamber design provides comfortable and effective treatment.

Features:

- Fully variable pressure range from 20 to 80 mmHg
- Lock button to prevent tampering
- Adjustable pause time
- Bilateral application when needed
- Four fully overlapping compression cells
- Color coding to eliminate error
- Quiet, safe operation
- Lightweight and portable
- **New!** Disposable garments for short term use



Pump Specifications

CYCLE	Sequential
PRESSURE RANGE	20-80 mmHg
SIZE	10.2x5.1x3.9 inches
WEIGHT	5.1 lb.
ELECTRICAL	115V, ~50/60Hz.

Quality you can count on

Manufacturing certified compliant with Good Manufacturing Requirements, ISO 13485:2003, and CAN/CSA ISO 13458:2003.

The Petite Basic System™

is approved by the FDA and UL.

