



Lympha Press USA

Lympha Press Optimal™

- Pretherapy™ based on the principles of manual lymph drainage
- Torso and abdominal therapy
- Bilateral or unilateral treatment
- 12-zone gradient sequential compression
- Therapy lock, for safe and easy use

The freedom of lymphedema treatment at home, with a therapist's touch.



Complete abdominal and torso therapy with our modular garments

Lympha Press Optimal™ is a home therapy system for lymphedema and venous insufficiency. It combines the gentle touch of manual lymph drainage with the proven efficacy of Lympha Press® gradient sequential compression, in a system that adapts to patient needs and is safe and easy to use.

Lympha Press Optimal™

Designed in accordance with the principles of Manual Lymph Drainage, Lympha Press Optimal™ is the most versatile, yet easy to use home care lymphedema treatment system available today.



Treatment modes:

- **Pretherapy™** – based on the principles of manual lymph drainage. This mode clears the proximal lymphatics prior to the main treatment cycle.
- **Sequential** – Lympha Press® quality 12-zone calibrated pneumatic compression therapy, clinically proven effective.
- **Peristaltic** – Wave-type calibrated gradient massage, for patients with focal discomfort.



MODES AND PRESSURE ARE LOCKED IN FOR SAFE PATIENT USE AT HOME.

Indications:

- Primary lymphedema
- Secondary lymphedema
- Venous insufficiency
- Venous stasis ulcers
- Dysfunction of the muscle pump

Specifications:

Model	Lympha Press Optimal™ 1201-AP
Pressure range	20 to 90 mmHg
Size	15×32×38.5 cm (6"×12 ½"×15¼")
Weight	6.5 kg (14.3 lbs)
Electrical	115V 50/60 Hz



FULL BILATERAL CAPABILITY:
Bilateral or unilateral therapy, with complete abdominal treatment.

Lympha Press USA Ltd. 232 Park Avenue, Manalapan, NJ 07726
Tel: 1-888-596-7421, 732-792-9677, Fax: 732-792-9745
info@lympha-press.com www.lympha-press.com