




Lympha Press USA



# LYMPHA PANTS™ II



UNIQUE modular design  
LEG SLEEVES zip on to torso garment  
TREATS one leg or two  
EASY to wear – just like a pair of pants!



## Lymphedema therapy with torso and abdominal treatment designed on the principles of manual lymph drainage

Leg lymphedema can be associated with lymphedema in the torso, groin, or abdomen. Until now, manual lymph drainage was the only option for torso therapy.

**Lympha Pants™ II - Torso treatment based on the principles of manual lymph drainage.** Designed in conjunction with leading manual lymph drainage therapists, **Lympha Pants II** treat the entire lower body, including the groin and abdomen up to the chest.

**Lympha Pants™ II** can be used to treat one leg or two simultaneously, with complete therapy of the entire abdomen and groin\*.

Use **Lympha Pants™ II** for:

- Leg lymphedema with abdominal, groin or torso edema
- History of genital lymphedema
- Bilateral leg lymphedema
- Unilateral leg lymphedema with spillover to the unaffected limb
- Leg lymphedema that has not responded to other forms of treatment



Use on one leg or two, together with complete abdominal therapy


\* We recommend Lympha Press leg sleeves for treatment of uncomplicated unilateral lymphedema.

# LYMPHA PANTS™ II

## Features:

- Complete abdominal, groin and leg therapy
- 12-24 overlapping chambers, to prevent reflux and promote directional flow
- Treats genital lymphedema and prevents its development in patients who may be susceptible
- Easy to put on and take off, just like a pair of jeans
- Easy to clean, for good health and skin care
- Non-stretch fabric, for correct compression every time
- Full bilateral capability

**Lympha Pants™ II are fully adjustable, for patients big and small.**



**Lympha Press Optimal™** - Calibrated gradient compression with Pretherapy™ based on the principles of manual lymph drainage. Pressure, therapy mode and treatment time are all adjustable to meet individual needs, with patient lock and easy on/off operation for home use. Medicare coded E0652.



**Lympha Press Mini®** - Calibrated gradient compression, with adjustable pressure in three zones, for comfortable and safe home therapy. Medicare coded E0652.

**Lympha Press USA Ltd.** 232 Park Avenue, Manalapan, NJ 07726  
Tel: 1-888-596-7421, 732-792-9677, Fax: 732-792-9745  
info@lymphapress.com www.lymphapress.com